

2D Electrophoresis

Trypsin in-gel Digest (PSC 537)

Contact: Hyuk-Kyu Seoh 404-413-5379; hseoh@gsu.edu

General Information There are numerous protocols for in-gel protein digestion - they are all variations on the same theme. This is the protocol that Kyu developed. The TFA has a limited shelf life. There is conflicting information about how long to use TFA once opened. Replace the TFA if it turns yellow or after 4 - 6 months. The Supplies section has order information about a 0.5 ml ampule of TFA, so you don't have to toss too much.

One of the biggest problems with protein yield after in-gel digestion is when the peptides stick to the plastic - tips, tubes or plates. Included in the Buffers and Supplies section are some low-bind products that might help. Also, be careful not to dry the peptides completely. This can't be avoided with a microplate, though. So, again, think about using a low-bind plate.

Protein digestion 5/1/07

Notes: • Perform reagent prep and protocol in hood using nitrile gloves • Thoroughly rinse all reagent tubes with ddH₂O

1. Process immediately after gel-plugs are picked.
2. Prepare first 3 reagent. Keep reconstituted trypsin at -20°C.
3. If plugs were collected using an Ettan Spotpicker, remove ddH₂O - be careful not to smash the plugs. If your gel was destained with 40% MeOH/10% Acetic acid, proceed to step 7.
4. Add 400µl 50% MeOH/10% acetic acid. Incubate 60 min and remove solution.
5. Repeat step 4.
6. Optional for storage: add 10 – 15 µl 50% MeOH/10% Acetic acid to store the gel-plug; store o/n or days at 4°C.
7. Remove all liquid and add 400µl 50% ACN/100 mM NH₄CO₃. Incubate 1 hour with agitation. Remove liquid.
8. Add 50µl 100% ACN. Incubate for 20 min. Remove liquid.
9. Dry plugs in speed vac (may need to remove rotor) for 1 h with no heat. Dried gel-plugs can be stored for days at -80°C)

10. When ready, make 1 X trypsin (20 ng/ μ l) from 10 X stock with 25 mM NH_4CO_3 (pH 8.0). Add small amount of 1 X trypsin solution to gel plug. Use 1 μ l to 7 μ l , depending on amount of protein in the gel-plug.
 - a. The amount of trypsin is usually 60 - 120 ng/gel plug depending on the amount of protein in the gel.
 - b. Keep the trypsin amount sub-molar ratio to protein.
 - c. Bring final volume to 30 - 40 μ l with 25 mM NH_4CO_3 (pH 8.0)
11. Seal plate and incubate at 37°C incubator overnight (14 – 16 hours).
12. Make up last two reagents.
13. Remove solution and put in new fresh plate/tube.
14. Add 30 - 40 μ l 50% acetonitrile/0.2% TFA. Incubate 30min and transfer to the same plate/tube as above.
15. Add 30 - 40 μ l 50% acetonitrile/0.1% TFA. Incubate 30min and transfer to the same plate/tube as above.
16. Dry solution in speedvac (no heat) until almost dry, around 1h or less. Check periodically to prevent drying to completion. If using tubes, try to NOT dry to completion. Peptide fragments tend to stick to plastic (use low bind tubes). Dry until about 10 μ l are left.
17. If you are going to clean up your sample with Zip-Tips, dry to 3 -5 μ l and resuspend to 10 μ l with 0.1% TFA.
18. Proceed to ziptip or freeze samples at -20C.

Buffers and Solutions

50% MeOH/ 10% acetic acid – 8ml

3.2ml ddH₂O

4ml MeOH

0.8ml acetic acid

100 mM ammonia bicarb / 50% Acetonitrile – 8ml

3.84ml ddH₂O

4ml ACN

800 μ l 1 M ammonium bicarb (pH 8.0)

200ug/ml trypsin/Promega reconstitution buffer = 10 X stock
to 100 µl reconstitution buffer (200ng/ µl)
mix and aliquot 5µl/tube

25 mM ammonia bicarb - 5 ml
125 µl 1M ammonium bicarb (pH 8.0)
4.875 ml ddH₂O

50% acetonitrile / 0.2% TFA – 4ml
2ml ddH₂O
2ml acetonitrile
8ul 100% TFA

50% acetonitrile / 0.1% TFA – 4ml
2ml ddH₂O
2ml acetonitrile
4ul 100% TFA

Supplies

1. Trypsin

Quantity/Unit: 100 ug/vial (5 X 20ug)

Vendor: Promega

Catalog#: V511A

2. TFA (Trifluoroacetic acid)

Quantity/Unit: 1 Pkg (5 X 0.5 ml)

Vendor: Aldrich

Catalog#: 308897-1PAK

3. Siliconized microfuge tubes, 0.5 ml (\$50.35/PK)

Quantity/Unit: 500/PK

Vendor: Applied Biosystems

Catalog#: AM12350

4. Siliconized microfuge tubes, 1.5 ml (\$40/PK)

Quantity/Unit: 250/PK

Vendor: Applied Biosystems

Catalog#: AM12450

5. Low Retention microfuge tubes, 0.6 ml (\$25/PK)

Quantity/Unit: 1000 (10 bags of 100)

Vendor: CLP, Continental Lab Products

Catalog#: 3435.S3

6. Low Retention microfuge tubes, 1.5 ml (\$14/PK)

Quantity/Unit: 500

Vendor: CLP, Continental Lab Products
Catalog#: 3445.S3

7. Low Retention pipet tips, Racked and Pre-sterilized, 200 ul(\$25.50/PK)

Quantity/Unit: 10 trays of 96

Vendor: CLP, Continental Lab Products

Catalog#: 2102.YS

8. Low Retention pipet tips, Racked and Pre-sterilized, 10 ul(\$25.50/PK)

Quantity/Unit: 10 trays of 96

Vendor: CLP, Continental Lab Products

Catalog#: 2142.S

9. Low Bind microfuge plates, conical, Pre-sterilized

Quantity/Unit: 100

Vendor: PGC

Catalog#: 81-6664-56