

SUMMER INSTITUTE CALENDAR 2022

SUN	MON	TUE	WED	THU	FRI	SAT
						July 02
						Early Arrival Airport Arrivals and Check-in
July 03	04	05	06	07	08	09
Early Arrival Airport Arrivals and Check-in	Airport Arrivals and Check-in 6:00pm: 4th of July Celebrations	9:30am-12pm: Campus tour, Panther ID & ISSS Check-in 12-2pm Lunch 2:00-6:00pm, Shuttle to local grocery store	9:30am-11:30am ISSS, OII, & Housing Orientation & Presentation 2:30-4:30pm: Welcome Reception and Buddy Meet & Greet Event	Classes begin! 9-11:20am: Morning course 11:20am-2:00pm: Lunch break 1:30-4:30pm: BIOL4905 INTRO - TRAINING	9-11:20am: Morning course 11:20am-2:00pm: Lunch break 1:30-4:30pm: BIOL4905 DNA PREPARATION	Free Day
10	11	12	13	14	15	16
12:00-4:00pm: The World Coca-Cola and Georgia Aquarium	9-11:20am: Morning course 11:20am-2:00pm: Lunch break 1:30-4:30pm: BIOL4905 PROTEOMICS I	9-11:20am: Morning course 11:20am-2:00pm: Lunch break 1:30-4:30pm: BIOL4905 PROTEOMICS II	9-11:20am: Morning course 11:20am-2:00pm: Lunch break 1:30-4:30pm: BIOL4905 PROTEOMICS III 6:00-10:00pm: Atlantic Station Shopping & Movie (Sign-up)	9-11:20am: Morning course 11:20am-2:00pm: Lunch break 1:30-4:30pm: BIOL4905 PROTEOMICS IV ?	9-11:20am: Morning course 11:20am-2:00pm: Lunch break 1:30-4:30pm: BIOL4905 RNA PREPARATION	6:00-9:00pm: Dinner in America (Sign-up)
17	18	19	20	21	22	23
Free Day	9-11:20am: Morning course 11:20am-2:00pm: Lunch break 1:30-4:30pm: BIOL4905 qPCR & AUTOMATION	9-11:20am: Morning course 11:20am-2:00pm: Lunch break 1:30-4:30pm: BIOL4905 DNA SEQUENCING	MINI BREAK	9-11:20am: Morning course CDC TRIP 1:30-4:30pm: BIOL4905 MICROSCOPY / AFM	9-11:20am: Morning course 11:20am-2:00pm: Lunch break 1:30 - 4:30pm: BIOL4905 NEXT GEN SEQ. 5:30-7:30pm: Meet & Greet BBQ event @ The Commons	9:00am - 6:00pm: Outlet Mall
24	25	26	27	28	29	30
Free Day	9-11:20am: Morning course 11:20am-2:00pm: Lunch break 1:30-4:30pm: BIOL4905 MICROARRAY I	9-11:20am: Morning course 11:20am-2:00pm: Lunch break 1:30-4:30pm: BIOL4905 NANOSTRING	9-11:20am: Morning course 11:20am-2:00pm: Lunch break Notebook Day 1:30-4:30pm: BIOL4905 NANOSTRING	Last day of classes 9-11:20am: Morning course 11:20am-2:00pm: Lunch break 1:30-4:30pm: BIOL4905 FLOW CYTOMETRY	FINALS	Free Day
31	August 01	02	03	04		
Free Day	Activity Day at the Recreation Center (Sign-up)	Free Day	9:30-11:00am: Georgia Capitol Tour (Sign-up) 2:00-4:00pm: Closing Reception	Departures (check-out at 12:00pm)		

Note: Students may arrive prior to the program date with an extra charge of \$35 per night. Earliest day to check-in to University Commons is July 2.

Legend:

Orange: Courses Blue: Lunch Break Red: Sign-up events

SUMMER INSTITUTE CALENDAR 2022

SUN	MON	TUE	WED	THU	FRI	SAT
						July 02
						Early Arrival Airport Arrivals and Check-in
July 03	04	05	06	07	08	09
Early Arrival Airport Arrivals and Check-in	Airport Arrivals and Check-in 6:00pm: 4th of July Celebrations	9:30am-12pm: Campus tour, Panther ID & ISSS Check-in 12-2pm Lunch 2:00-6:00pm, Shuttle to local grocery store	9:30am-11:30am ISSS, OII, & Housing Orientation & Presentation 2:30-4:30pm: Welcome Reception and Buddy Meet & Greet Event	Classes begin! 9-11:20am: Morning course 11:20am-2:00pm: Lunch break 1:30-4:30pm: BIOL4905 INTRO - TRAINING	9-11:20am: Morning course 11:20am-2:00pm: Lunch break 1:30-4:30pm: BIOL4905 DNA PREPARATION	Free Day
10	11	12	13	14	15	16
12:00-4:00pm: The World Coca-Cola and Georgia Aquarium	9-11:20am: Morning course 11:20am-2:00pm: Lunch break 1:30-4:30pm: BIOL4905 PROTEOMICS I	9-11:20am: Morning course 11:20am-2:00pm: Lunch break 1:30-4:30pm: BIOL4905 PROTEOMICS II	9-11:20am: Morning course 11:20am-2:00pm: Lunch break 1:30-4:30pm: BIOL4905 PROTEOMICS III 6:00-10:00pm: Atlantic Station Shopping & Movie (Sign-up)	9-11:20am: Morning course 11:20am-2:00pm: Lunch break 1:30-4:30pm: BIOL4905 PROTEOMICS IV ?	9-11:20am: Morning course 11:20am-2:00pm: Lunch break 1:30-4:30pm: BIOL4905 RNA PREPARATION	6:00-9:00pm: Dinner in America (Sign-up)
17	18	19	20	21	22	23
Free Day	9-11:20am: Morning course 11:20am-2:00pm: Lunch break 1:30-4:30pm: BIOL4905 qPCR & AUTOMATION	9-11:20am: Morning course 11:20am-2:00pm: Lunch break 1:30-4:30pm: BIOL4905 DNA SEQUENCING	MINI BREAK	9-11:20am: Morning course CDC TRIP 1:30-4:30pm: BIOL4905 MICROSCOPY / AFM	9-11:20am: Morning course 11:20am-2:00pm: Lunch break 1:30 - 4:30pm: BIOL4905 NEXT GEN SEQ. 5:30-7:30pm: Meet & Greet BBQ event @ The Commons	9:00am - 6:00pm: Outlet Mall
24	25	26	27	28	29	30
Free Day	9-11:20am: Morning course 11:20am-2:00pm: Lunch break 1:30-4:30pm: BIOL4905 MICROARRAY I	9-11:20am: Morning course 11:20am-2:00pm: Lunch break 1:30-4:30pm: BIOL4905 NANOSTRING	9-11:20am: Morning course 11:20am-2:00pm: Lunch break 1:30-4:30pm: BIOL4905 NANOSTRING Notebook Day	Last day of classes 9-11:20am: Morning course 11:20am-2:00pm: Lunch break 1:30-4:30pm: BIOL4905 FLOW CYTOMETRY	FINALS	Free Day
31	August 01	02	03	04		
Free Day	Activity Day at the Recreation Center (Sign-up)	Free Day	9:30-11:00am: Georgia Capitol Tour (Sign-up) 2:00-4:00pm: Closing Reception	Departures (check-out at 12:00pm)		

Note: Students may arrive prior to the program date with an extra charge of \$35 per night. Earliest day to check-in to University Commons is July 2.

Legend:

Orange: Courses Blue: Lunch Break Red: Sign-up events

SUMMER INSTITUTE CALENDAR 2022


SUN	MON	TUE	WED	THU	FRI	SAT
						July 02
						Early Arrival Airport Arrivals and Check-in
July 03	04	05	06	07	08	09
Early Arrival Airport Arrivals and Check-in	Airport Arrivals and Check-in 6:00pm: 4th of July Celebrations	9:30am-12pm: Campus tour, Panther ID & ISSS Check-in 12-2pm Lunch 2:00-6:00pm, Shuttle to local grocery store	9:30am-11:30am ISSS, OII, & Housing Orientation & Presentation 2:30-4:30pm: Welcome Reception and Buddy Meet & Greet Event	Classes begin! 9-11:20am: Morning course 11:20am-2:00pm: Lunch break 1:30-4:30pm: BIOL4905 INTRO - TRAINING	9-11:20am: Morning course 11:20am-2:00pm: Lunch break 1:30-4:30pm: BIOL4905 DNA PREPARATION	Free Day
10	11	12	13	14	15	16
12:00-4:00pm: The World Coca-Cola and Georgia Aquarium	9-11:20am: Morning course 11:20am-2:00pm: Lunch break 1:30-4:30pm: BIOL4905 PROTEOMICS I	9-11:20am: Morning course 11:20am-2:00pm: Lunch break 1:30-4:30pm: BIOL4905 PROTEOMICS II	9-11:20am: Morning course 11:20am-2:00pm: Lunch break 1:30-4:30pm: BIOL4905 PROTEOMICS III 6:00-10:00pm: Atlantic Station Shopping & Movie (Sign-up)	9-11:20am: Morning course 11:20am-2:00pm: Lunch break 1:30-4:30pm: BIOL4905 PROTEOMICS IV ?	9-11:20am: Morning course 11:20am-2:00pm: Lunch break 1:30-4:30pm: BIOL4905 RNA PREPARATION	6:00-9:00pm: Dinner in America (Sign-up)
17	18	19	20	21	22	23
Free Day	9-11:20am: Morning course 11:20am-2:00pm: Lunch break 1:30-4:30pm: BIOL4905 qPCR & AUTOMATION	9-11:20am: Morning course 11:20am-2:00pm: Lunch break 1:30-4:30pm: BIOL4905 DNA SEQUENCING	MINI BREAK	9-11:20am: Morning course CDC TRIP 1:30-4:30pm: BIOL4905 MICROSCOPY / AFM	9-11:20am: Morning course 11:20am-2:00pm: Lunch break 1:30 - 4:30pm: BIOL4905 NEXT GEN SEQ. 5:30-7:30pm: Meet & Greet BBQ event @ The Commons	9:00am - 6:00pm: Outlet Mall
24	25	26	27	28	29	30
Free Day	9-11:20am: Morning course 11:20am-2:00pm: Lunch break 1:30-4:30pm: BIOL4905 MICROARRAY I	9-11:20am: Morning course 11:20am-2:00pm: Lunch break 1:30-4:30pm: BIOL4905 NANOSTRING	9-11:20am: Morning course 11:20am-2:00pm: Lunch break 1:30-4:30pm: BIOL4905 NANOSTRING Notebook Day	Last day of classes 9-11:20am: Morning course 11:20am-2:00pm: Lunch break 1:30-4:30pm: BIOL4905 FLOW CYTOMETRY	FINALS	Free Day
31	August 01	02	03	04		
Free Day	Activity Day at the Recreation Center (Sign-up)	Free Day	9:30-11:00am: Georgia Capitol Tour (Sign-up) 2:00-4:00pm: Closing Reception	Departures (check-out at 12:00pm)		

Note: Students may arrive prior to the program date with an extra charge of \$35 per night. Earliest day to check-in to University Commons is July 2.

Legend:

Orange: Courses Blue: Lunch Break Red: Sign-up events

Questions




Is this a good qPCR plot of data?


Explain your answer in your notebook.

Questions


EcoR I



Sma I



Pst I



Which of these Type II REN can you find in YOUR sequence (next page). If you find the site, after which base number in the sequence does it cut?

Explain your answer by providing your analysis in your notebook.



GSU Biology Core Facility




Supporting Life Sciences at GSU

http://biotech.gsu.edu/core_facility/index.html





Flow Cytometry

Profiling Cells






Cellular Functions

Mass Spectrometry



Atomic Force Microscopy


Imaging at the Ångström level




DNA Sequence Analysis:

Profiling DNA

Sanger Sequencing –
>800 base pairs/run



High Throughput Genomic Sequencing –
100,000,000,000 base pairs/run



RNA Expression

Microarray: Analysis

Profiling mRNA



Colour of pin-point dots demonstrates the presence / absence of gene sequences

Protein Expression

Proteomics

Profiling Proteins

2D Protein gel
Protein separation using
Electric charge and molecular weight






GSU Biology Core Facility

Supporting Life Sciences at GSU

http://biotech.gsu.edu/core_facility/index.html





Flow Cytometry Profiling Cells





Cellular Functions

Mass Spectrometry




Atomic Force Microscopy Imaging at the Ångström level




Sanger Sequencing –
>800 base pairs/run

DNA Sequence Analysis: Profiling DNA



High Throughput Genomic Sequencing –
100,000,000,000 base pairs/run



Microarray:Analysis Profiling mRNA



Colour of pin-point dots demonstrates the presence / absence of gene sequences



Proteomics

Profiling Proteins

2D Protein gel
Protein separation using
Electric charge and molecular weight



GSU Biology Core Facility



Supporting Life Sciences at GSU

http://biotech.gsu.edu/core_facility/index.html




Flow Cytometry

Profiling Cells





See effects of different drugs on Cell cycle



Apoptosis -programmed cell death

Cellular Functions

Mass Spectrometry



Proteomics




Profiling Proteins

2D Protein gel
Protein separation using Electric charge and molecular weight

Atomic Force Microscopy

Imaging at the Ångström level


Protein structure analysis




Sanger Sequencing –
>800 base pairs/run

DNA Sequence Analysis:

Profiling DNA



High Throughput Genomic Sequencing –
100,000,000,000 base pairs/run



Microarray: Analysis

Profiling mRNA



Colour of pin-point dots demonstrates the presence / absence of gene sequences



GSU Biology Core Facility




Supporting Life Sciences at GSU

http://biotech.gsu.edu/core_facility/index.html





Flow Cytometry

Profiling Cells





Cellular Functions

Mass Spectrometry



Atomic Force Microscopy

Imaging at the Ångström level



Protein Expression

Proteomics

Profiling Proteins

2D Protein gel
Protein separation using Electric charge and molecular weight


DNA

Replication


Sanger Sequencing –
>800 base pairs/run

DNA Sequence Analysis:

Profiling DNA



High Throughput Genomic Sequencing –
100,000,000,000 base pairs/run



RNA Expression

Microarray: Analysis

Profiling mRNA



Colour of pin-point dots demonstrates the presence / absence of gene sequences

An Overview of GeneChip[®] Technology -



John Houghton, PhD
GSU ABCore Facilities

What is Microarray?

Microarrays

circa 1991

(Schena et al. (1995) *Science* 270:467-70)

Probe DNA is attached to solid support
plastic beads, glass slide, nylon
or chip

RNA is labeled (usually indirectly)

Arrays can detect

mRNA

microRNA

Methylation

SNP

High throughput


10,000s of specific probes

Measure global gene expression,
SNP calls, LOH, amplification, methylation etc




GeneChip® vs. Spotted Arrays


- Affymetrix GeneChip® Arrays use oligonucleotides
- Oligos are built on a solid support



- Spotted arrays utilize nucleic acids made in solution
- Solutions are then “spotted” onto a solid support
- Competitive Hybridization




About Affymetrix




Historically....

DNA spotting

- DNA spotting usually uses multiple pins
- DNA in microtiter plate
- DNA usually PCR amplified
- Oligonucleotides can also be spotted




Afymetrix cornered the market




Sequence Database



Research Tools




Expression Variability



Sequence Variability

About Affymetrix



GeneChip Human Genome U133A 2.0 Array

The GeneChip® Human Genome U133A 2.0 Array is a single array representing 14,500 well-characterized human genes that can be used to explore human biology and disease processes.

Part #	Description	Unit Size	Your Price (USD)	Qty
900471	Human Genome U133A 2.0 Array	contains 2 arrays	Please Inquire	
900468	Human Genome U133A 2.0 Array	contains 6 arrays	Please Inquire	
900469	Human Genome U133A 2.0 Array	contains 30 arrays	Please Inquire	

- Product Description
- Technical Documentation
- Required/Related Products

The GeneChip® Human Genome U133A 2.0 Array is a single array representing 14,500 well-characterized human genes that can be used to explore human biology and disease processes. New design and reduced feature size mean that you can use smaller sample volumes than the previous HG@U133A Array without compromising performance.

- Provides coverage of well-substantiated genes in the transcribed human genome on a single array
- Analyzes the expression level of 18,400 transcripts and variants, including 14,500 well-characterized human genes
- Comprised of more than 22,000 probe sets and 500,000 distinct oligonucleotide features
- Use the Power of the Probe Set and get multiple independent measurements for each transcript that deliver the greatest accuracy and reproducibility of any microarray platform
- All probe sets represented on the GeneChip® Human Genome U133A Array are identically replicated on the GeneChip Human Genome U133A 2.0 Array

Array Profile

Sequences used in the design of the array were selected from GenBank®, dbEST, and RefSeq. The sequence clusters were created from the UniGene database (Build 133, April 20, 2001) and then were refined by analysis and comparison with a number of other publicly available databases including the Washington University EST trace repository and the University of California, Santa Cruz Golden-Path human genome database (April 2001 release).

Instrument and Software Requirements

- GeneChip® Scanner 3000, enabled for High-Resolution Scanning*
- GeneChip® Command Console® Software (AGCC) including the GeneChip® Scanner 3000 High-Resolution Scanning Patch

*GeneChip Scanner 3000 High-Resolution Update is standard on all instruments shipped starting in September 2003 with serial number series 502. Previous versions, serial number series 501, will require the 00-0110 GeneChip Scanner 3000 High-Resolution Update to be installed.

For more information, please review the [data sheet](#) (pdf, 169 KB).

For research use only. Not for use in diagnostic procedures.

Expression Arrays 2006

Feline Gene 1.0 ST Array	30 arrays	Please Inquire	
Guinea Pig Gene 1.0 ST Array	Contains 6 arrays	Please Inquire	
Guinea Pig Gene 1.0 ST Array	Contains 30 arrays	Please Inquire	
Marmoset Gene 1.0 ST Array	6 arrays	Please Inquire	
Marmoset Gene 1.0 ST Array	30 arrays	Please Inquire	
Medicago Gene 1.0 ST Array	6 arrays	Please Inquire	
Medicago Gene 1.0 ST Array	30 arrays	Please Inquire	
Ovine Gene 1.0 ST Array	6 arrays	Please Inquire	
Ovine Gene 1.0 ST Array	30 arrays	Please Inquire	
Porcine Gene 1.0 ST Array	6 arrays	Please Inquire	
Porcine Gene 1.0 ST Array	30 arrays	Please Inquire	
Rabbit Gene 1.0 ST Array	Contains 6 arrays	Please Inquire	
Rabbit Gene 1.0 ST Array	30 arrays	Please Inquire	
Rhesus Gene 1.0 ST Array	6 arrays	Please Inquire	
Rhesus Gene 1.0 ST Array	30 arrays	Please Inquire	
Rice (Cn) Gene 1.0 ST Array	30 arrays	Please Inquire	
Rice (Cn) Gene 1.0 ST Array	6 arrays	Please Inquire	
901997	Rice (Jp) Gene 1.0 ST Array	6 arrays	Please Inquire
901996	Rice (Jp) Gene 1.0 ST Array	30 arrays	Please Inquire
901991	Rice (US) Gene 1.0 ST Array	30 arrays	Please Inquire
901992	Rice (US) Gene 1.0 ST Array	6 arrays	Please Inquire
902002	Soybean Gene 1.0 ST Array	6 arrays	Please Inquire
902001	Soybean Gene 1.0 ST Array	30 arrays	Please Inquire
902299	Tomato Gene 1.0 ST Array	6 arrays	Please Inquire
902300	Tomato Gene 1.0 ST Array	30 arrays	Please Inquire
901956	Zebra Finch Gene 1.0 ST Array	30 arrays	Please Inquire
901957	Zebra Finch Gene 1.0 ST Array	6 arrays	Please Inquire
902007	Zebrafish Gene 1.0 ST Array	6 arrays	Please Inquire
902006	Zebrafish Gene 1.0 ST Array	30 arrays	Please Inquire

Rat Human

Murine

Arabidopsis

C. elegans

Drosophila

Usa

Yeast

Eukarya

Enabling the Genetic Revolution

Food testing




Livestock diagnostics or grading




Environmental testing



Agricultural biotech



Basic Research



Identity testing



Individualized medicine



Human diagnostics

Understanding Information -Human Genome Project


Gene Functions at a Basic Level

- Gene Identification
 - Which genes are important and in which tissues?
- Pathway Characterization
 - Define relationships between genes
- Regulation
 - Examine motifs on a global scale

Specific Applications in Healthcare & Pharma



- Tumour Typing
 - Use expression patterns to complement classical histology to identify classes of tumors and predict disease development
- Drug Response
 - Monitor impact of a therapeutic on disease state or toxicological effect

Understanding Cancer




Distinguishing Between Leukaemias

Pediatric **A**cute **M**yeloid **L**eukaemia (**AML**)
T or B cell **A**cute **L**ymphoblastic **L**eukaemia (**ALL**)
Mixed **L**ineage **L**eukamias (**MLL**)



ALL



AML



MLL

Which class of
leukemia?

Individual Variations in Drug Metabolism

Cytochrome P450, Detoxification Enzymes

Cytochrome P450 (CYP) enzymes are a superfamily of **mono-oxygenases** that are found in all kingdoms of life, and which show extraordinary diversity in their reaction chemistry. In mammals, these enzymes are found primarily in the membranes of the **endoplasmic reticulum (microsomes)** within **liver cells** (hepatocytes), as well as many other cell types. These enzymes use **haeme iron** to **oxidise molecules**, often making them more water-soluble for clearance.

They achieve this by either adding or unmasking a polar group. In general, the reaction catalysed by these enzymes can be summarised as:



Individual Variations in Drug Metabolism

“Intestinal cytochrome P450 proteins play an important role in the biotransformation of drugs and may significantly limit their oral absorption.”

Drug Metabolism and Disposition, June 2008 vol. 36 no. 6 1039-1045

[Int J Clin Pharmacol Res. 2003;23\(1\):31-5.](#)

Genetic polymorphism of cytochrome P450 enzymes in Asian populations: focus on CYP2D6.

[Kitada M¹.](#)

⊕ Author information

Abstract

Published studies demonstrate that significant ethnic differences can exist in the metabolism of some drugs. These differences are caused by cytochrome P450 polymorphisms and result in the potential for wide interpatient and interethnic variability in adverse events. One of the most common of these cytochrome P450 polymorphisms is related to the CYP2D6 isozyme. Many classes of commonly used drugs are metabolized by CYP2D6, creating the potential for significant adverse events. Due to the variety of genetic polymorphisms among Asian populations, this article focuses on this group rather than on other ethnic populations and discusses the clinical importance of genetic polymorphisms with regard to potential drug interactions. Polymorphism of CYP2D6 can either increase the rate of drug elimination (ultrametabolizers, leading to faster metabolic clearance potentially resulting in reduced effectiveness and need for higher doses) or decrease drug metabolism (poor metabolizers, which may increase the potential for drug interactions and adverse events). Although the CYP2D6 poor metabolizer phenotype is less frequent in Asian than in Western populations (e.g. about 1% in Thai, Chinese and Japanese populations and up to 4.8% in Indians versus 5-10% in Caucasians), the increased prevalence of the CYP2D6*10 allele in Asians does have an impact on drugs metabolized by CYP2D6. Enzyme activity is reduced, potentially increasing circulating drug doses and increasing the risk for drug interactions. Thus, in Asian populations it may be important to optimize pharmacotherapy either by assessing patients' CYP2D6 genotype, or by prescribing medications that are not metabolized by this isozyme.

PMID: [14621071](#)

[PubMed - indexed for MEDLINE]

Individual Variations in Drug Metabolism



CYP450 genes metabolize more than 90% of commercially available drugs



- Poor Metabolizers**
- Two mutant alleles
 - No enzyme activity



- Intermediate Metabolizers**
- One reduced activity allele
 - One null allele



- Extensive Metabolizers**
- At least one normal allele



- Ultrarapid Metabolizers**
- Multiple functional alleles
 - Excess enzymatic activity

Individual Variations in Drug Metabolism



CYP450 genes metabolize more than 90% of commercially available drugs



- Poor Metabolizers**
- Two mutant alleles
 - No enzyme activity



- Intermediate Metabolizers**
- One reduced activity allele
 - One null allele



- Extensive Metabolizers**
- At least one normal allele



- Ultrarapid Metabolizers**
- Multiple functional alleles
 - Excess enzymatic activity

Individual Variations in Drug Metabolism



CYP450 genes metabolize more than 90% of commercially available drugs



- Poor Metabolizers**
- Two mutant alleles
 - No enzyme activity



- Intermediate Metabolizers**
- One reduced activity allele
 - One null allele



- Extensive Metabolizers**
- At least one normal allele



- Ultrarapid Metabolizers**
- Multiple functional alleles
 - Excess enzymatic activity

Individual Variations in Drug Metabolism



CYP450 genes metabolize more than 90% of commercially available drugs



- Poor Metabolizers**
- Two mutant alleles
 - No enzyme activity



- Intermediate Metabolizers**
- One reduced activity allele
 - One null allele




- Extensive Metabolizers**
- At least one normal allele



- Ultrarapid Metabolizers**
- Multiple functional alleles
 - Excess enzymatic activity

Individual Variations in Drug Metabolism



CYP450 genes metabolize more than 90% of commercially available drugs

What is GeneChip® Technology?



GeneChip® System




**Probe
Arrays
(chips)**



**Fluidics
Station**




Scanner




**Software
Data Analysis**

GeneChip[®] Technology?




49 Chips per Wafer



400 Chips per Wafer

GeneChip Probe Arrays



5"
5"
49 - 400
Chips/wafer

1.28cm

Image of Hybridized Probe Array


GeneChip[®] Array Advantages

- Assume fixed array size, 1.28 x 1.28 cm

<u>Feature Size</u>	<u>Features/Chip</u>	<u>Genes/Chip</u>
100 μm	16,384	409
50 μm	65,538	1,638
24 μm	284,444	7,111
20 μm	409,600	12,800
18 μm	506,944	~ 22,500
10 μm	1,600,000	>200,000

GeneChip® Expression Analysis


Hybridization and Staining




GeneChip[®] Expression Analysis

Hybridization and Staining


Probe Array



Hybridized Array




Labeled cRNA Target



Staining



Streptavidin-
phycoerythrin
conjugate



+



=





Avidin

Biotin conjugate

Avidin-Biotin Complex


1.28 cm
1.28 cm
Actual size of
GeneChip[®] array




Millions of DNA strands built up in each location

500,000 locations on each GeneChip[®] array


Actual strand = 25 base pairs




**Sample RNA fragments (purple)
washed over DNA probe array (green)**




**Sample RNA fragments (purple)
hybridized to DNA probe array (green)**




**Sample RNA fragments (purple)
washed over DNA probe array (green)**




**C does not stick to another C,
so no match is made**



We know there was no match because there is no fluorescent RNA bound to the probe.



**Sample RNA fragments (purple)
hybridized to DNA probe array (green)**




***Labeled RNA (purple) binds to the matching probes in the single feature.**


***The feature will now fluoresce**

GeneChip® vs. Spotted Arrays


- Affymetrix GeneChip® Arrays use oligonucleotides
 - Oligos are built on a solid support




- Spotted arrays utilize nucleic acids made in solution
 - Solutions are then “spotted” onto a solid support
 - Competitive Hybridization




Proven Manufacturing Systems



Wafer Prep



Photolithography




Cartridge Assembly



**Finished GeneChip
Probe Array**


Proven Manufacturing Systems



Wafer Prep



Cartridge Assembly




Photolithography




**Finished GeneChip
Probe Array**

Proven Manufacturing Systems



Wafer Prep



Photolithography




Cartridge Assembly




**Finished GeneChip
Probe Array**

Proven Manufacturing Systems



Wafer Prep



Photolithography




Cartridge Assembly



**Finished GeneChip
Probe Array**


Photolithographic Synthesis




Photolithography




Photolithographic Synthesis




Photolithographic Synthesis




How Probes Are Synthesized




How Probes Are Synthesized




How Probes Are Synthesized



How Probes Are Synthesized



How Probes Are Synthesized



How Probes Are Synthesized




C A T A T
A G C T G
T T C C G
~~~~~

# Combinatorial Synthesis

*Any N-mer can be synthesized in  $4 \times N$  steps*


## Example:

Sixteen 20-mers



< 80-step  
synthesis

$\sim 10^{12}$  20-mers



< 80-step  
synthesis

# Applications of Gene Expression

- Basic Research
  - Cell Cycle
  - Signaling Pathways
  - Regulation of cell differentiation
- 
- Genetic Basis of Disease
  - Genetic changes in cancer; classification
  - Metabolic diseases
  - Aging-related biological pathologies
  - Immune system pathologies
  - Infectious diseases
- 
- Target Discovery and Drug Development
  - Analyze disease models to discover drug targets
  - Profile drug candidates using expression data
  - Analyze drug toxicity in various model systems


# GeneChip<sup>®</sup> Array Advantages

- Multiple Indicators for the Same Target Ensures:

- Specificity
- Quantitative accuracy
- Low false positive rate
- High sensitivity


# GeneChip<sup>®</sup> Expression Analysis

## Hybridization and Staining



# GeneChip<sup>®</sup> Array Advantages

## Specificity




Gene



Oligo probes for Increased Gene Specificity


# GeneChip<sup>®</sup> Array Advantages

## Specificity



# GeneChip<sup>®</sup> Array Advantages

## Specificity





# GeneChip<sup>®</sup> Array Advantages

## Specificity

Gene 1



— Oligo probes for Gene 1

Gene 2



— Oligo probes for Gene 2

cDNA  
for Gene 1



cDNA  
for Gene 2




Spotted  
cDNAs

# Detecting Change in Gene Expression

*Yeast grown in different conditions*

Rich Medium



Minimal Medium



Genes

a = RPL2A

b = TIP1

c = BAP2

d = VAP1

e = YBR147W

f = SUP46


Wodicka, L., et al. 1997  
*Nature Biotechnology*, 15: 1359-1367



# GeneChip<sup>®</sup> Array Advantages


## Specific Hybridization

Called  
“Present”




## Non-specific/Cross Hybridization

Called  
“Absent”





# 25-mer PM-MM Probe Pairs Offer Increased Specificity



Discrimination between target and stray signal at low (<8pM) target concentrations facilitated by the use of MM probes.

# Positioning of primer probes can be important



Images of Hybridized Probe  
Arrays

Over 500,000 different probes  
complementary to genetic  
information of interest


# Sensitivity vs. Specificity

- **Sensitivity**
  - Identifying low abundance transcripts
  - Tolerate some miscalls to achieve greater sensitivity
  - Avoid false negatives
- **Specificity**
  - Accuracy of detection
  - Tolerate missing some calls to achieve accuracy
  - Avoid false positives

# Reproducibility


Independent cell growth and prep

Sample 1



Signal Intensities: 1 13,400  
2 1,280  
3 1 (absent)  
4 1,840  
5 1,700


Sample 2



Signal Intensities: 1 11,670  
2 1,250  
3 9 (absent)  
4 2,010  
5 1,450




# Reproducibility



# Reproducibility


Same sample, different arrays

Sample 1




**Signal Intensities:** 13,400  
1,280  
1 (absent)  
1,840  
1,700

Sample 1



**Signal Intensities:** 13,090  
1,250  
10 (absent)  
1,750  
1,430


# Reproducibility




# Differences Across Commercial Platforms

Increased Expression

Decreased Expression





# Differences Across Commercial Platforms






# Differences Across Commercial Platforms



# Sources of Variability




| Experiment    | Label     | Cell                 | Fluor  | Channel | ExpType |
|---------------|-----------|----------------------|--------|---------|---------|
| T. wild mouse | H2K-K26C2 | SGC-26005.108        | 144.90 | 1       | 144.90  |
| T. wild mouse | H2K-K207  | N4T-1027.07.983      | 144.90 | 2       | 144.90  |
| T. wild mouse | H2K-K209  | N4T-1029.05.609      | 144.90 | 3       | 144.90  |
| T. wild mouse | H2K-K284  | SGC-26042.0.193      | 144.90 | 3       | 144.90  |
| T. wild mouse | H2K-K219  | SGC-26049.0.182      | 144.90 | 3       | 144.90  |
| T. wild mouse | H2K-K205  | N4T-1027.04.021      | 144.90 | 3       | 144.90  |
| T. wild mouse | H2K-K202  | T18B-48642.0.201.948 | 144.90 | 3       | 144.90  |
| T. wild mouse | H2K-K288  | SGC-26052.0.102      | 144.90 | 3       | 144.90  |
| T. wild mouse | H2K-K276  | T18B-48642.0.193.005 | 144.90 | 3       | 144.90  |
| T. wild mouse | H2K-K285  | SGC-26051.0.102      | 144.90 | 3       | 144.90  |
| T. wild mouse | H2K-K282  | SGC-26049.0.187      | 144.90 | 3       | 144.90  |
| T. wild mouse | H2K-K208  | N4T-1027.04.021      | 144.90 | 3       | 144.90  |
| T. wild mouse | H2K-K277  | SGC-26043.0.171      | 144.90 | 3       | 144.90  |
| T. wild mouse | H2K-K242  | N4T-1027.01.463      | 144.90 | 3       | 144.90  |
| T. wild mouse | H2K-K286  | SGC-26043.0.176      | 144.90 | 3       | 144.90  |
| T. wild mouse | H2K-K271  | SGC-26050.0.144      | 144.90 | 3       | 144.90  |
| T. wild mouse | H2K-K203  | N4T-1029.05.733      | 144.90 | 3       | 144.90  |
| T. wild mouse | H2K-K289  | SGC-26050.0.180      | 144.90 | 3       | 144.90  |
| T. wild mouse | H2K-K283  | SGC-26050.0.181      | 144.90 | 3       | 144.90  |
| T. wild mouse | H2K-K207  | N4T-1027.02.803      | 144.90 | 3       | 144.90  |
| T. wild mouse | H2K-K206  | N4T-1029.05.734      | 144.90 | 3       | 144.90  |
| T. wild mouse | H2K-K289  | SGC-26051.0.171      | 144.90 | 3       | 144.90  |
| T. wild mouse | H2K-K273  | T18B-48642.0.176     | 144.90 | 3       | 144.90  |
| T. wild mouse | H2K-K204  | N4T-1029.05.102      | 144.90 | 3       | 144.90  |
| T. wild mouse | H2K-K201  | T18B-48642.0.189     | 144.90 | 3       | 144.90  |
| T. wild mouse | H2K-K200  | N4T-1027.01.772      | 144.90 | 3       | 144.90  |
| T. wild mouse | H2K-K204  | SGC-26049.0.108      | 144.90 | 3       | 144.90  |

- **Biology**
  - The main source of variability
- **Sample preparation**
  - Technical variability depends on method and operator
- **Probe array analysis**
  - Standardized; relatively little variability
- **Data analysis**



# GeneChip<sup>®</sup> System Work Flow



# GeneChip Expression Analysis



Absolute Analysis



- Detection (qualitative)
- Signal (quantitative)




Comparison Analysis



- Change (qualitative)
- Signal Log Ratio (quant.)

# Genotyping or Gene Expression Monitoring



# GeneChip<sup>®</sup> Mapping 10K Probe Design



probes = 25 bases

Perfect Match

Mismatch

Perfect Match

Mismatch



Allele 'A'



Allele 'G'

# SNP Tiling Strategy

SNP 0 Position

  
G / A

TAGCCATCGGTA **N** GTA **C** TCAATGATCAGCT

---

|              |          |                       |           |          |
|--------------|----------|-----------------------|-----------|----------|
| ATCGGTAGCCAT | <b>C</b> | CAT <b>G</b> AGTTACTA | PM Allele | <b>A</b> |
| ATCGGTAGCCAT | <b>G</b> | CAT <b>G</b> AGTTACTA | MM Allele | <b>A</b> |
| ATCGGTAGCCAT | <b>T</b> | CAT <b>G</b> AGTTACTA | PM Allele | <b>B</b> |
| ATCGGTAGCCAT | <b>A</b> | CAT <b>G</b> AGTTACTA | MM Allele | <b>B</b> |

---

# SNP Tiling Strategy

## Single nucleotide polymorphisms: aging and diseases

B Bessenyei <sup>1</sup>, M Márka, L Urbán, M Zeher, I Semsei

Affiliations + expand

PMID: 15547317 DOI: [10.1007/s10522-004-2567-y](https://doi.org/10.1007/s10522-004-2567-y)

### Abstract

Differences of more than 3 million nucleotides can be seen comparing the genomes of two individuals as a result of single nucleotide polymorphism (SNP). More and more SNPs can be identified and it seems that these alterations are behind of several biological phenomena. Personal differences in these nucleotides result for example in elevated disease susceptibilities, that is, certain nucleotides are more frequent in patients suffering from different diseases comparing to the healthy population. SNPs may cause substantial alterations in the cells, e.g. the enzyme activity of the respective gene changes, but in other cases the effects of the SNPs are not so pronounced. Later results indicate that SNPs can be rendered to individuals living a longer life than the average. Perhaps these results will not directly lead to the lengthening of the maximal life span; however, genes that play an important role in the aging process could be identified. In this respect SNPs are important factors in determining the information level of the cells of individuals which determines the maximal life span (I. Semsei On the nature of aging. Mech. Ageing Dev . 2000; 117: 93-108), in turn SNP is one of the factors that determine the aging process. Since there are certain age-related diseases, the discovery and the description of the SNPs as a function of age and diseases may result in a better understanding of the common roots of aging and those diseases.

# SNP Tiling Strategy

## Heart disease

**Heart disease**, in SNPedia as well as for the entry in [Wikipedia](#), is a catch-all term including medical classifications such as coronary artery disease, myocardial infarction, atherosclerosis, etc. Heart disease overall is the #1 cause of death in developed countries, typically accounting for up to 40% of all deaths.

Many SNPs have been associated with increased risk for one or more types of heart disease. Before listing many of them, though, it is worth emphasizing that the risks associated with these SNPs add "relatively little to the current capacity of traditional, non-genetic risk factors to identify individuals with a high propensity to develop heart disease"<sup>[1]</sup>. This is generally true of most SNPs associated with other diseases as well.

SNPs and genes associated with altered risk for heart disease include the following:

- The most highly replicated associations to **heart disease** have been to SNPs in the chromosome 9p21 region. SNPs in this region include:
  - [rs2383206](#)
  - [rs10757278](#)
  - [rs2383207](#)
  - [rs10757274](#)

# SNP Tiling Strategy

## Heart disease

---

**Heart disease**, in SNPedia as well as for the entry in [Wikipedia](#), is a catch-all term including medical classifications such as coronary artery disease, myocardial infarction, atherosclerosis, etc. Heart disease overall is the #1 cause of death in developed countries, typically accounting for up to 40% of all deaths.

Many SNPs have been associated with increased risk for one or more types of heart disease. Before listing many of them, though, it is worth emphasizing that the risks associated with these SNPs add "relatively little to the current capacity of traditional, non-genetic risk factors to identify individuals with a high propensity to develop heart disease"<sup>[1]</sup>. This is generally true of most SNPs associated with other diseases as well.

SNPs and genes associated with altered risk for heart disease include the following:

- The most highly replicated associations to **heart disease** have been to SNPs in the chromosome 9p21 region. SNPs in this region include:
  - [rs2383206](#)
  - [rs10757278](#)
  - [rs2383207](#)
  - [rs10757274](#)



# SNP Tiling Strategy

## Heart disease

**Heart disease**, in SNpedia as well as for the entry in [Wikipedia](#), is a catch-all term including medical classifications such as coronary artery disease, myocardial infarction, atherosclerosis, etc. Heart disease overall is the #1 cause of death in developed countries, typically accounting for up to 40% of all deaths.

Many SNPs have been associated with increased risk for one or more types of heart disease. Before listing many of them, though, it is worth emphasizing that the risks associated with these SNPs add "relatively little to the current capacity of traditional, non-genetic risk factors to identify individuals with a high propensity to develop heart disease"[1]. This is generally true of most SNPs associated with other diseases as well.

SNPs and genes associated with altered risk for heart disease include the following:

- The most highly replicated associations to **heart disease** have been to SNPs in the chromosome 9p21 region. SNPs in this region include:
  - [rs2383206](#)
  - [rs10757278](#)
  - [rs2383207](#)
  - [rs10757274](#)

## rs2383206

[rs10757274](#) and [rs2383206](#) can significantly increase the risk of [heart disease](#)[1]. About one in every four Caucasians are thought to carry the variants, and their risk of coronary heart disease is increased by 30 to 40%. [rs10757278](#) in the same region has been linked to [diabetes](#) [2]. The chromosomal region where these SNPs are located is 9p21, and has no known genes.

a [blog post](#) about investigating [rs10757274](#) and [rs2383206](#)

[[PMID 18048766](#)] This SNP was also associated with increased risk for [coronary artery disease](#) in a Korean population.

[[PMID 18066490](#)] Also found to be significant in a study of 416 Italian myocardial infarction patients.

A study of 1,000+ patients with early-onset angiographic [coronary artery disease](#) (CAD) concluded that [rs2383206](#)(G) was associated with an adjusted odds ratio of 1.39 (CI: 1.05-1.85) for (A;G) heterozygotes and 1.73 (CI: 1.26-2.37) for (G;G) homozygotes. This SNP alone accounted for 21% of the [population attributable fraction](#) and was independent of traditional risk factors, myocardial infarction risk, and the extent of disease.[[19033013?dopt=Abstract PMID 19033013](#)]

[[PMID 19559344](#)] Genetic variants on chromosome 9p21 and ischemic stroke in Chinese

|                 |                 |                                       |   |
|-----------------|-----------------|---------------------------------------|---|
| Orientation     | plus            |                                       |   |
| Stabilized      | plus            |                                       |   |
| Geno            | Mag             | Summary                               | ↕ |
| (A;A)           | 1               | normal                                |   |
| (A;G)           | 2               | 1.4x increased risk for heart disease |   |
| (G;G)           | 3               | 1.7x increased risk for heart disease |   |
| Reference       | GRCh38 38.1/141 |                                       |   |
| Chromosome      | 9               |                                       |   |
| Position        | 22115027        |                                       |   |
| Gene            | CDKN2B-AS1      |                                       |   |
| is a            | snp             |                                       |   |
| is              | mentioned by    |                                       |   |
| dbSNP           | rs2383206       |                                       |   |
| dbSNP (classic) | rs2383206       |                                       |   |
| dbSNP (classic) | rs2383206       |                                       |   |

# SNP Tiling Strategy

## Heart disease

**Heart disease**, in SNpedia as well as for the entry in [Wikipedia](#), is a catch-all term including medical classifications such as coronary artery disease, myocardial infarction, atherosclerosis, etc. Heart disease overall is the #1 cause of death in developed countries, typically accounting for up to 40% of all deaths.

Many SNPs have been associated with increased risk for one or more types of heart disease. Before listing many of them, though, it is worth emphasizing that the risks associated with these SNPs add "relatively little to the current capacity of traditional, non-genetic risk factors to identify individuals with a high propensity to develop heart disease"[1]. This is generally true of most SNPs associated with other diseases as well.

SNPs and genes associated with altered risk for heart disease include the following:

- The most highly replicated associations to **heart disease** have been to SNPs in the chromosome 9p21 region. SNPs in this region include:
  - [rs2383206](#)
  - [rs10757278](#)
  - [rs2383207](#)
  - [rs10757274](#)

## rs2383206

[rs10757274](#) and [rs2383206](#) can significantly increase the risk of [heart disease](#)[1]. About one in every four Caucasians are thought to carry the variants, and their risk of coronary heart disease is increased by 30 to 40%. [rs10757278](#) in the same region has been linked to [diabetes](#) [2]. The chromosomal region where these SNPs are located is 9p21, and has no known genes.

a [blog post](#) about investigating [rs10757274](#) and [rs2383206](#)

[PMID 18048766] This SNP was also associated with increased risk for [coronary artery disease](#) in a Korean population.

[PMID 18066490] Also found to be significant in a study of 416 Italian myocardial infarction patients.

A study of 1,000+ patients with early-onset angiographic [coronary artery disease](#) (CAD) concluded that [rs2383206](#)(G) was associated with an adjusted odds ratio of 1.39 (CI: 1.05-1.85) for (A;G) heterozygotes and 1.73 (CI: 1.26-2.37) for (G;G) homozygotes. This SNP alone accounted for 21% of the [population attributable fraction](#) and was independent of traditional risk factors, myocardial infarction risk, and the extent of disease.[[19033013?dopt=Abstract PMID 19033013](#)]

[PMID 19559344] Genetic variants on chromosome 9p21 and ischemic stroke in Chinese

|             |      |                                       |  |
|-------------|------|---------------------------------------|--|
| Orientation | plus |                                       |  |
| Stabilized  | plus |                                       |  |
| Geno        | Mag  | Summary                               |  |
| (A;A)       | 1    | normal                                |  |
| (A;G)       | 2    | 1.4x increased risk for heart disease |  |
| (G;G)       | 3    | 1.7x increased risk for heart disease |  |

Reference [GRCh38 38.1/141](#)

Chromosome 9

Position 22115027

Gene [CDKN2B-AS1](#)

is a [snp](#)

is [mentioned by](#)



dbSNP [rs2383206](#)

dbSNP [rs2383206](#)

(classic)

dbSNP [rs2383206](#)

# GeneChip Mapping 10K Assay




- A single-primer,
- A single-amplification, and
- A cost-effective assay...
- ...for simultaneously **Genotyping**

# GeneChip<sup>®</sup> Mapping 10K Hybridization Patterns

**AA**

**BB**

**AB**



PHASE : INTERPRETATION  
TWO

SEIDMAN *with* Leizer

